



WORKSTATION SPECIFICATIONS

SYNAPSE 

FUJIFILM
Value from Innovation

Learn how **Synapse® 5, the fastest PACS in the medical industry**, delivers more through flexible architecture and diverse workstation specifications.

Client Specifications

Version 4.4.3 Client

Operating System	Vendor	Browser	CPU	RAM	Comments
Windows 7 Ultimate with SP1 32- and 64-bit	Microsoft	Internet Explorer 8, 9, 10, or 11	Clinical: Dual Core (minimum) Diagnostics: Quad Core (minimum)	Clinical: 8 GB (minimum) Diagnostics: 16 GB (minimum)	Standard VGA and PCIe video card-compliant
Windows 7 Professional with SP1 32- and 64-bit	Microsoft	Internet Explorer 8 ,9 ,10, or 11	Clinical: Dual Core (minimum) Diagnostics: Quad Core (minimum)	Clinical: 8 GB (minimum) Diagnostics: 16 GB (minimum)	Standard VGA and PCIe video card-compliant
Windows 7 Home Premium with SP1 32- and 64-bit	Microsoft	Internet Explorer 8 ,9 ,10, or 11	Clinical: Dual Core (minimum) Diagnostics: Quad Core (minimum)	Clinical: 8 GB (minimum) Diagnostics: 16 GB (minimum)	Standard VGA and PCIe video card-compliant Home Premium edition does not allow domain login. IE should be configured to prompt for password.
Windows 8.1 32- and 64-bit	Microsoft	Internet Explorer 11	Clinical: Dual Core (minimum) Diagnostics: Quad Core (minimum)	Clinical: 8 GB (minimum) Diagnostics: 16 GB (minimum)	Standard VGA and PCIe video card-compliant

Client Specifications

Version 5.2 Client

Operating System	Vendor	Browser	CPU	RAM	Comments
Windows 7 Ultimate with SP1 64-bit	Microsoft	Internet Explorer 11 Google Chrome Browser	Clinical: Dual Core (minimum) Diagnostics: Quad Core (minimum)	Clinical: 8 GB (minimum) Diagnostics: 16 GB (minimum)	Standard VGA and PCIe video card-compliant
Windows 7 Professional with SP1 64-bit	Microsoft	Internet Explorer 11 Google Chrome Browser	Clinical: Dual Core (minimum) Diagnostics: Quad Core (minimum)	Clinical: 8 GB (minimum) Diagnostics: 16 GB (minimum)	Standard VGA and PCIe video card-compliant
Windows 7 Home Premium with SP1 64-bit	Microsoft	Internet Explorer 11 Google Chrome Browser	Clinical: Dual Core (minimum) Diagnostics: Quad Core (minimum)	Clinical: 8 GB (minimum) Diagnostics: 16 GB (minimum)	Standard VGA and PCIe video card-compliant Home Premium edition does not allow domain login. IE should be configured to prompt for password.
Windows 8.1 32- and 64-bit	Microsoft	Internet Explorer 11 Google Chrome Browser	Clinical: Dual Core (minimum) Diagnostics: Quad Core (minimum)	Clinical: 8 GB (minimum) Diagnostics: 16 GB (minimum)	Standard VGA and PCIe video card-compliant
Windows 10 32- and 64-bit	Microsoft	Internet Explorer 11 Google Chrome Browser	Clinical: Dual Core (minimum) Diagnostics: Quad Core (minimum)	Clinical: 8 GB (minimum) Diagnostics: 16 GB (minimum)	Standard VGA and PCIe video card-compliant



More Speed
Less Bandwidth
Tighter Security
Comprehensive Enterprise Imaging
Improved User Efficiency

Instantly access even the largest datasets with the most popular browsers, using less bandwidth throughout the enterprise.

Learn how **Synapse 5** can improve performance in your organization.

If you would like to schedule a product demonstration or have a representative contact you, please click below.

CONTACT US